



- ✓ Better Data
- ✓ Better Decisions
- ✓ Better Performance




The Equivital™ TnR Product Range

July 2011 - Sept 2011

- ✔ Better Data
- ✔ Better Decisions
- ✔ Better Performance

The Equivital™
TnR Product Range



Equivital™ TnR Products combine world class human performance monitoring capability with robust clinical grade data quality.

Designed for use by real people in real environments, the products keep usability and flexibility at the fore, meeting broad end-user monitoring requirements.

Equivital™ technologies have pioneered mobile human monitoring for almost a decade. The TnR Product range brings these proven Equivital™ technologies to users in the form of new sensor systems and software viewing & analysis tools.

Equivital™ EQ02 LifeMonitor

The Equivital™ LifeMonitor is a multi-parameter body worn sensor capable of logging and transmitting physiological, psychophysiological and biomechanical data.



EQ02 Sensor Belt

Product code: EQ-02-B1-1-S(SIZE), EQ-02-B2-1-S(SIZE)*

The EQ02 Sensor belt features a system of sensors and electrodes embedded in novel textiles. It connects with the Sensor Electronics Module (SEM) to provide mobile monitoring capability.

Designed to be comfortable when worn for extended periods, the EQ02 Sensor belt enables high quality measurement of physiology even when the wearer is engaged in strenuous physical activity and rapid movement.

Available sizes:

1-9 (lower chest circumference: 29 – 47 in / 74 – 119 cm).

Reference:

Item 1.1 within price list.



EQ02 SEM

Product code: EQ-02-SEM-007

The EQ02 SEM is the worlds leading multi parameter ambulatory monitoring device that senses, records, monitors and intelligently processes data captured from the human and is able to transmit this over a wireless or wired interface.

The EQ02 SEM weights 38gm, is IP67 certified when connected with a sensor belt and has class 1 Bluetooth interface.

Reference:

Item 2 within price list.



EQ02 M-Dock

Product code: EQ-ACC-MD-1

Enables simultaneous charging and 2-way data transfer communication with up to six SEM's.

Five M-Docks can be daisy chained to support 30 concurrent connections via USB.

Reference:

Item 3 within price list.

*This model of the EQ02 Sensor belt features an external connector for wired ancillary sensors or enables connection to an additional external power pack.

- ✓ Better Data
- ✓ Better Decisions
- ✓ Better Performance

Software Platforms & Communications

Data can be viewed in real time or retrospectively using Equivital™ viewing and analysis software tools.

The TnR Software range includes Equivital manager, a data management and configuration tool, and the eqView Software suite for web, local and mobile access to Equivital sensor data.



Product Code: EQ-SW-12

A PC application used for configuration of SEM's and data management.

Reference: Software Item 5.5



Product Code: EQ-SW-EQV-P-1

A local or networked PC application used for monitoring subjects, viewing live data, trends and replaying logged data.

Reference: Software Item 5



Product Code: EQ-SW-EQV-EXP-1

A web service, which receives live data from a subject using eqView Mobile, allowing the user to monitor subjects remotely.

Reference: Software Item 5.2






Product Code: EQ-SW-EQV-MOB-A-1

A mobile application used for viewing live data from a subject locally. eqView Mobile can also send data via GPRS, SMS or WI-FI to eqView Explorer. (Currently versions for RIM and Android enabled phones available).


Reference: Software Item 5.3

Equivital™ eqView Professional Configurations



- ✓ Better Data
- ✓ Better Decisions
- ✓ Better Performance

Equivital™ eqView Explorer based System Configuration



Equivital™ LifeMonitor Accessories



EQ02 Single Lead

Product Code: EQ-ACC-SL-1

Enables simultaneous charging and 2-way data transfer with a single SEM.

Reference: Accessories 3.1



EQ-GSR Galvanic Skin Response Sensor

Product Code: EQ-ACC-034

Wrist mounted sensor for measurement of galvanic skin response. (Suitable for use with EQ02 Belt +).

Reference: Accessories 3.5



EQ Core Temperature Capsule & Dermal Patch Activator

Product Code: EQ-ACC-029

Simply activates the core temperature capsule and the dermal temperature patch. (Requires the use of a radio receiver to test capsule and patch activation).

Reference: Accessories 3.8



Nonin Xpod

Product Code: EQ-ACC-042

Low power connection to SEM for wired oxygen saturation probes and various compatible Nonin sensors*.

*Other sensors are not supplied by Equivital™.

Reference: Accessories 3.3



Dermal Patch

Product Code: EQ-ACC-048

Self adhesive dermal temperature patch measures skin temperature and transmits data to the SEM in real time.

Reference: Accessories 3.6



Nonin Ipod

Product Code: EQ-ACC-041

Oxygen saturation probe with SEM connector and finger clip.

Reference: Accessories 3.4



Core Temperature Capsule

Product Code: EQ-ACC-023

Ingestible core temperature capsule measures core body temperature and transmits readings to the SEM in real time.

Reference: Accessories 3.7



Equivital™ Hub

Product Code: EQ-ACC-32

The Equivital™ Hub is a Bluetooth access point enabling communication with up to 18 SEMs from a PC or LAN in real time. Networkable via WiFi or Ethernet connection. Operates with eqView Professional software.

Reference: Communications 4

Further information

If you require any further information about Hidalgo or Equivital™ please contact us at the address below:

Hidalgo Limited

Unit F, Trinity Court,
Buckingway Business Park
Swavesey
Cambridge
CB24 4UQ
United Kingdom

T +44(0) 1954 233430

F +44(0) 1954 233431

E info@equivital.co.uk

W www.equivital.co.uk

Copyright © Hidalgo Limited. All rights reserved. All other company and product names maybe trademarks of their respective companies. While every effort is made to make sure the information shown is accurate, Hidalgo does not accept any liability for any errors or mistakes that may arise. Specifications and other information may be subject to change without notice. All performance figures and other data contained in this document may vary by application.

Equivital™ by Hidalgo - a Jaltex Group Company

HIDA3330-TNRPR-02-1.1