

Ultra Mirage™ II Nasal Mask

COMPONENTS CARD

Ultra Mirage II Mask System—ROW

(complete system with cushion and headgear)

Standard	16545
Large	16546
Shallow	16547
Shallow-wide	16578

Ultra Mirage II Mask System—UK

(complete system with cushion and headgear)

Standard	16321
Large	16322
Shallow	16323
Shallow-wide	16324

Ultra Mirage II Frame System

(no cushion and headgear)

Frame	16595
-------	--------------

Ultra Mirage II Mask System—Australia

(complete system with cushion and headgear)

Standard	16579
Large	16580
Shallow	16581
Shallow-wide	16582

Headgear
(not shown)
16118 Small
16733 Standard
16119 Large

Forehead pads
60123 (pack)
60124 (10 pack)

Forehead support with pads
16590

Mask frame
16560

Elbow clip
16572
16573 (10 pack)

Swivel
16565

Cushion clip
16563

Cushion
16556 (standard)
16557 (large)
16558 (shallow)
16735 (shallow-wide)

Vent cover
16575

Headgear clip
16569 (2 pack)
16734 (10 pack)

Ports cap
16570 (2 pack)
16571 (10 pack)

Elbow assembly
includes: Elbow clip, vent
cover, elbow and swivel
16566

Masque Ultra Mirage™ II

ECLATÉ DU MASQUE

Masques Ultra Mirage

Masque Ultra Mirage - standard	16545	Ensemble rigide	16595
Masque Ultra Mirage - large	16546	(comprend le support frontal	
Masque Ultra Mirage - shallow	16547	& l'entourage rigide;	
Masque Ultra Mirage - shallow/wide	16578	sans bulle, ni harnais)	

Harnais
(non représenté)
16118 Petit
16733 Standard
16119 Large

Clip de la bulle
16563

Bulle
16556 (standard)
16557 (large)
16558 (shallow)
16735 (shallow/wide)

Clip d'attache
du harnais
16569 (par 2)
16734 (par 10)

Bouchons d'entrée
16570 (par 2)
16571 (par 10)

Tampon frontal
16561 (par 2)
16574 (par 10)

Support frontal
16590

Entourage rigide
16560

Clip du coude
16572
16573 (par 10)

Pièce
Pivotante
16565

Couvercle orifices
de ventilation
16575

Coude rotatif
(Comprend le clip du coude, le
couvercle orifices de ventilation,
le coude et la pièce pivotante)
16566

Ultra Mirage™ II Maske

KOMPONENTENKARTE

Ultra Mirage Maske (Komplettmaske)

Standard	16545
Large	16546
Shallow	16547
Shallow-wide	16578

Maskenhaltersystem

(ohne Maskenkissen und Kopfband)

Maskenhalter	16595
--------------	--------------

Kopfband
(nicht abgebildet)
16118 Small
16733 Standard
16119 Large

Stützring
16563

Maskenkissen
16556 (standard)
16557 (large)
16558 (shallow)
16735 (shallow-wide)

Kopfbandclip
16569 (2 Stck.)
16734 (10 Stck.)

Verschlussstöpsel
16570 (2 Stck.)
16571 (10 Stck.)

Stirnpolster
16561 (2 Stck.)
16574 (10 Stck.)

Stirnbauteil
16590

Maskenhalter
16560

Sprengring
16572
16573 (10 Stck.)

Drehadapter
16565

Luftauslasskappe
16575

Kniestück (Beinhaltet:
Sprengring, Luftauslasskappe,
Kniestück und Drehadapter)
16566

